YGC Agar (Yeast-Glucose-Chloramphenicol)

Catalog #	Description
3555489	YGC Agar (Yeast-Glucose-Chloramphenicol) , ready-to-use, 100 ml x 6 bottles
3564104	YGC Agar (Yeast-Glucose-Chloramphenicol) , dehydrated, 500 g

For laboratory use only.

Intended Use

Medium used for the enumeration of yeast and mold in the analysis of food products, especially in dairy products.

Principle

The nutrient substances provided by yeast extract and the use of glucose as an energy source promote the growth of yeast and mold. The presence of chloramphenicol, a thermostable broad-spectrum antibiotic, inhibits the growth of non-target microorganisms.

Theoretical Composition

Base Medium

Yeast extract	5 g
Glucose	20 g
Chloramphenicol	0.1 g
Agar	16 g
Distilled water	1,000 ml
Final pH at 25°C	= 6.6 ± 0.2

Shelf Life and Storage

Store ready-to-use media at 2–25°C. Store dehydrated medium at 15–25°C in carefully sealed bottles in a cool, dry place.

Required Materials Not Supplied

This list is not exhaustive.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer
- Water-bath capable of maintaining ±1°C

Supplies

- 100 ml Pyrex bottles with autoclave proof stoppers
- Sterile Petri dishes (Ø = 90 mm)
- Sterile pipettes (0.1 ml, 1 ml, etc)
- Sterile Pasteur pipettes
- Sterile spreaders

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Avoid prolonged heating during melting

- The medium may look frothy after solidification in bottles. It nevertheless keeps all its qualities even if its appearance changes slightly during melting and shaking
- The time lapse between the end of preparation of the stock solution (or the 10⁻¹ dilution in the case of a solid product) and the moment when the dilutions come into contact with the culture medium must not exceed 15 min
- Petri dishes should be handled with care in order to avoid the dispersion of mold spores
- In cases where heavy contamination by Gram negative bacteria is suspected (particularly in meat or raw fish), it is recommended that a solution of Gentamicin be added to the medium (final concentration = 100 mg/L), previously cooled to 44–47°C
- For SDS product safety information and certificate of analysis, visit bio-rad.com

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

Dehydrated Medium Preparation

- Shake bottle before use
- Dissolve 41.1 g of powder in 1 L of sterile distilled water
- Wait 5 min, then mix until a homogenous suspension is obtained
- Heat gently, agitating frequently, then bring to a boil until completely dissolved
- Dispense into 100 bottles and sterilize in an autoclave at 121 ± 3°C for 15 min

Reconstitution ratio: 41.1 g/L (500 g of powder makes 12.1 L of medium)

Sample Preparation

- Prepare sample according to the standard method applicable to the product concerned

Inoculation and incubation

Surface inoculation

- Deposit 0.1 ml of sample and/or its decimal dilutions onto the surface of the dried agar
- Spread the inoculum quickly but carefully using a sterile spreader and allow to dry

Depth inoculation:

- Transfer 1 ml of each sample and/or dilution to be analyzed in a sterile Petri dish
- Pour 15 ml of tempered (44–47°C) agar into the dish. Homogenize by swirling and allow to solidify on a cool, level surface

Incubation:

- Invert plates and incubate at 25 ± 1°C for 5 days, depending on the standard

Reading and Interpretation

- Enumerate the colonies on each plate and express the results according to ISO 7218

Note: If necessary, distinguish colonies of yeast and mold from colonies of bacteria by morphological characteristics under microscopic examination

References

ISO 6611:2004 [IDF 94:2004]. Milk and milk products — Enumeration of colony-forming units of yeast and/or moulds — Colony-count technique at 25 degree C.

ISO 7218:2007/AMD1:2013. Microbiology of food and animal feeding stuffs — General requirements and guidance for microbiological examinations.

Mossel, DAA., Visser, M, and Mengerink, WHJ (1962). A comparison of media for the enumeration of molds and yeasts in food and beverages. *Practical Microbiology* 11:650-653.

NF V08-059 November 2002. Microbiology of food animal feeding stuffs — Enumeration of yeasts and moulds by colony-count technique at 25 °C — Routine method.

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Revision History

Release date	Document number	Change
July 2022	5098 Ver A	<ul style="list-style-type: none">- Major change- New document design- Document number change — previous version: V4 12-08-11

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

YGC Agar (Yeast-Glucose-Chloramphenicol)

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